


BeHear[®] NOW Personalizable Hearing Amplifier




Sound enhancement technologies by:

Worldwide Hearing Facts

~15% HEARING LOSS




<5% USE HEARING AIDS




Good News for the Hearing Impaired



Popular Hearing Complaints





What can I hear when I have:

Normal Hearing

Mild Hearing Loss


Moderate Hearing Loss

Mod/Severe Hearing Loss

Severe Hearing Loss

Profound Hearing Loss

Reasons for Low Hearing Aid Uptake




Performance




DSP sound enhancement technologies

Since
2002






Powering >50M
products

www.alango.com



Simple, Immediate, Self-Managed



Neck Loop Form Factor Houses High Quality Components



1. 2x2 low-noise microphones
2. 13mm HiFi speakers
3. Powerful 120MHz DSP
4. Long life battery



Built-in Hearing Assessment



Personalization at Your Fingertips



“Best hearing point” – **user** finds it!

Bonus Technologies: ListenThrough™



Stay safe and
alert when
listening to
streamed audio

Before



After



Slow down
speech during
mobile calls

Bonus Technologies: HearLink™



BeHear

Now



HearLink
W&H transmitter



Any standard television set
with analog audio out port

**HearLink is sold separately*

Addressing the Full Hearing Spectrum



Thank You

Sound Enhancement by Alango Technologies, Ltd.

info@wearandhear.com

www.wearandhear.com