


BeHear[®] ACCESS

Personalizable Hearing Amplifier with Telecoil




Worldwide Hearing Facts

~15% HEARING LOSS




<5% USE HEARING AIDS




Good News for the Hearing Impaired



Popular Hearing Complaints





What can I hear when I have:

Normal Hearing

Mild Hearing Loss


Moderate Hearing Loss

Mod/Severe Hearing Loss

Severe Hearing Loss

Profound Hearing Loss

Reasons for Low Hearing Aid Uptake




Performance




DSP sound enhancement technologies

Since
2002




Powering >50M
products

www.alango.com



Simple, Immediate, Quality Hearing

Wear&Hear
Life Sounds Better




Neck Loop Form Factor Houses High Quality Components


1. 2x2 low-noise microphones
2. 2 telecoil receivers (in ACCESS)
3. 13mm HiFi speakers
4. Powerful 120MHz DSP
5. Long life battery




Dedicated Smartphone App



Pair & Personalize



Select Hearing Preset






Customize in Real-time

Built-in Hearing Assessment




Personalization at Your Fingertips



“Best hearing point” – **user** finds it!

Built-in Telecoil Support

- 2 orthogonally-positioned telecoil receivers deliver optimal sound
- Exclusive and Transparent modes
- ADA-compliant



Bonus Technologies: ListenThrough™



Stay safe and
alert when
listening to
streamed audio

Before



After



Slow down
speech during
mobile calls

Bonus Technologies: HearLink®



Low Latency
<40ms



Supports
2 Headsets



High Quality
Stereo Sound

**HearLink is sold separately*

Bonus Technologies: HearLink® PLUS



Low Latency
<40ms



Long Range
Wireless



Supports
2 Headsets



High Quality
Stereo Sound

**HearLink PLUS is sold separately*

Addressing the Full Hearing Spectrum



Mobile Calls & Audio



Public Places



Live Music



TV Listening



Live Conversations



BeHear NOW vs BeHear ACCESS



Feature	BeHear NOW	BeHear ACCESS
Personalizable hearing enhancement	+	+
Real-time customization (W&H BeHear App)	+	+
EasyListen™ & ListenThrough™ functionality	+	+
HearLink® & HearLink® PLUS transmitter support	+	+
Highly fashionable headset with low-profile buttons	+	×
Telecoil receivers for use with hearing loops	×	+
Large, tactile control buttons on headset	×	+
Magnetized charging cradle	×	+
Hearing mode battery life (hours)	12+	13+
Price (MSRP)	\$249	\$349

Thank You

Sound Enhancement by Alango Technologies, Ltd.

info@wearandhear.com

www.wearandhear.com

